



Standard ACT Features

- Patient-friendly functionality for improved compliance
- Up to 30 days of retrievable ECG Memory
- Remotely programmable
- Automatically recorded and transmitted daily recording

Unique Features

- 3 channel ECG for more specific and sensitive Diagnostic Yields
- Provision of 24-48 hour Holter with retrievable ECG memory
- Ability for QT and ST segment analysis
- 6 Hour Sensor Memory

Best Use

Symptomatic/Asymptomatic AF
 Syncope
 Conductive Disorders
 Tachycardia, Bradycardia
 Non-compliance
 Medication Titration
 Pediatrics
 Cardiac Rehabilitation

General Description

The ACT III Platinum offers broad functionality with its high-performance multi-channel Ambulatory Cardiac Telemetry in one simple-to-use system. The ACT III Platinum system offers up to 30-days of real-time ECG monitoring.

With the provision of QT analysis and ST deviation that can measure ischemic changes, the ACT III Platinum will also provide more sensitive and specific data for initial or early detection of arrhythmia in patients that have minimal or atypical symptoms.

The ACT III Platinum utilizes a sophisticated algorithm and advanced wireless and web-based technologies for enhanced data collection, analysis, transmission and reporting.

- Automatically records and sends episodes of:
 - Atrial fibrillation (any rate)
 - Bradycardia
 - Tachycardia
 - Pause
- Manual activation button



LifeWatch ACT Services

We offer the most advanced technologies, comprehensive monitoring services, and clinical expertise to deliver the highest diagnostic yield.

